Digital Scale
with  Bluetooth[®]

Model : HD-351BT

(Eng) Instruction Manual

Read this Instruction Manual carefully and keep it for future reference.

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
1. Precaution

Precautions


To recalibrate the scale, gently tap the Personal Key with your toe (Fig.1). Do not place your hand on the scale while recalibrating (Fig.2); if you do "-0.2," or "-0.4," will be displayed. Using your toe, tap the Personal Key again to display "0.0"

Care and Maintenance

- Place the scale in an area free from heating equipment, high humidity, or extreme temperature changes.
- Never submerge the scale in water or spill chemical liquids on it. To clean, use a dampened cloth and mild detergent.
- Avoid excessive impact or vibration to the unit.

2. Part Names


- ❶ Measuring Platform
- ❷ LCD Display
- ❸ Battery Compartment on Bottom
- ❹ Unit Mode Switch on Bottom
- ❺ Personal Key
- ❻ Bluetooth OUTPUT Key

3. Feature

- Output weight data using Bluetooth technology.
- Previous weight using personal key.
- Digital accuracy and reliability.
- Classic design fits any décor.
- Durable construction; non-skid surface.
- "Lo" message appears when batteries need replacing.

4. Preparation Before Use

⚠ Inserting the Batteries

Open the battery cover on the back of the measuring platform. Insert the supplied AA-batteries as indicated. Be sure to use only alkaline batteries.

Note : Be sure that the polarity of the batteries is set properly. If the batteries are incorrectly positioned the fluid may leak and damage floors. If you do not intend to use this unit for a long period of time, it is advisable to remove the batteries before storage.

Please note that the included batteries from factory may have decreased energy levels over time.

⚠ Positioning the Monitor

Place the measuring platform on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

Note : To avoid possible injury, do not step on the edge of the platform.


5. Switch The Weight Mode

Note: For use only with models containing the weight mode switch

The units can be changed using the switch on the back of the scale as shown in the figure.

lb (pounds)

kg (kilograms)

st-lb (stone pounds)

6. Operation without Bluetooth

Instructions for Use


Tap the personal key with your toe.

The display will show personal number (Fig.3).

Wait for the display to show "0.0" (Fig.4), then step on the scale (Fig.5).

Current reading in upper display and previous reading in lower display.

- Note :**
- OL will appear on the display when the scales weight capacity is exceeded.
 - "0.0" disappears if you do not step on the scale within 60 seconds.
 - The unit will shut off if you tap the OUTPUT key, except --- while measuring.


7. Operation with Bluetooth

Initial Pairing

1-1) Before activating the scale, make sure the corresponding Bluetooth device is turned on.

1-2) Press "OUTPUT" key to turn on the scale, as noted by a single beep.


1-3) Within 30 seconds, activate the communication sequence within the Bluetooth device to initiate the pairing process with the scale. If 30 seconds is exceeded an error message will appear and the process must be repeated (press the "OUTPUT" key).


1-4) The pairing process is completed on the corresponding Bluetooth device and the scale, as noted by a double beep.

1-5) The display will change to the zero calibration display.

1-6) The display will then change to the waiting display. Step onto the scale and measure weight.

1-7) The weight result will flash after measurement is completed. The weight data will then be sent to the corresponding Bluetooth device.

1-8) To continue with more weightings, press "OUTPUT" key (return to step 1-5).

1-9) The scale will automatically shut off when ---

*The user remains on the scale for 30 seconds.

*The user get off from the scale, and then after 30 seconds.

After Initial Pairing


2-1) Before activating the scale, make sure the corresponding Bluetooth device is turned on.

2-2) Press "OUTPUT" key to turn on the scale, as noted by a single beep.

The scale will find the previous corresponding Bluetooth device.


* If unable to establish pairing or if no response from the corresponding Bluetooth device within 30 seconds, an error will show on the scale (step 3-1).


2-3) A successful pairing operation between scale and the corresponding Bluetooth device will be noted by a double beep.

2-4) The display will change to the zero calibration display.

2-5) The display will then changed to the waiting display. Step onto the scale and measure weight.


2-6) The weight result will flash after measurement is completed. The weight data will then be sent to the corresponding Bluetooth device.

2-7) To continue with more weightings, press "OUTPUT" key (Step 2-4).


2-8) The scale will automatically shut off when ---

* The user remains on the scale for 30 seconds.

* The user get off from the scale, and then after 30 seconds.

**Error Process**


3-1) Failure to "PAIR" display.

3-2) Display will change to zero calibration display.

3-3) The display will then change to the waiting display. Step onto the scale and measure weight, or to reattempt Bluetooth pairing press "OUTPUT" key (step 1-2).

3-4) The weight result will flash after measurement is completed.

* The weight result data will not be sent to the corresponding Bluetooth device.

* To reattempt Bluetooth pairing, press "OUTPUT" key (step 1-2).

3-5) The scale will automatically shut off when ---

* The user remains on the scale for 30 seconds.

* The user get off from the scale, and then after 30 seconds.

8. Output Data Format**Specification**

Signal Speed	9600 bps
Start Bit	1 bit
Data Bit Length	8 bits
Parity	None
Stop bit	1 bit
PIN	0000

Output Data

	Heade	Output data (ASCII code)	Substance
Control Data	{0	Fix to 16	2bytes fixed
Control Data	~0	Fix to 1	1bytes fixed
Model No.	MO	"HD-351BT"	10bytes fixed
Weight Data	kg	Wk	3 - 5bytes variable length
	lb	Wp	3 - 5bytes variable length
	St-lb	Ws	5 - 7bytes variable length X:st Y:lb
Checksum	CS	XX	2bytes fixed

Note : • The data are divided with commas (,) for each data.

• The terminator (end of data) is CR (ASCII code ^(zero)0D), LF (ASCII code ^(zero)0A).

[Output example (in the case to lb unit)]

Control Data	Control Data	Model No	Weight Data	Checksum	CR	LF															
{0	,	16	,	~0	,	1	,	MO	,	"HD-351BT"	,	Wp	,	xxx.x	,	CS	,	xx	,	^(zero) 0D	^(zero) 0A

9. Specifications

Maximum Capacity	440 lb / 200 kg / 31st 6lb
Graduation	0.2 lb / 0.1 kg / 0.2 lb
Power Supply	DC 6V (LR6 - AA Alkaline battery x 4 included)
Power Consumption	60 mA maximum
During Transmission	120 mA or less
When On	60 mA or less
When Display is Off	50 µA or less
Temperature Range for Use	41 to 95°F (5 to 35°C), 80% or less humidity (non condensation)
Temperature Range for Storage	41 to 95°F (5 to 35°C), 80% or less humidity (non condensation)

10. Bluetooth Specifications

Product Reference	Bluetooth Module
MITSUMI Part Number	WML-C46ASR
Application Specification	Bluetooth Ver. 2.0+EDR
Frequency Range	2.402 – 2.480 GHz
Transmission Scheme	FESS
Modulation Method	GFSK
Output Power	Class 2
Communications Range	Approx. 33 ft / Approx. 10 m (differs by usage environment)

 **Bluetooth**® The Bluetooth word mark and logo are owned by the Bluetooth SIG, Inc. and any use of such marks by Tanita Corporation is under license. Any other trademarks or trade names mentioned are the property of their respective owners.

"This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation." This device contains

FCC-ID POOWML-C46

Modifications

The FCC requires the user to be notified that any changes or modifications made to this device that are not expressly approved by Tanita Corporation may void the user's authority to operate the equipment.

This device has been approved for use in USA.

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TANITA Corporation
TANITA Corporation of America, inc.
Illinois 60005 U.S.A.