
TANITA®
Monitoring Your Health

BF-522W

Body Fat Monitor / Scale

Tanita UK Ltd.

The Barn, Philpots Close, Yiewsley,
Middlesex, UB7 7RY UK
Tel: +44 (0)1895 438577
Fax: +44 (0)1895 438511
www.tanita.co.uk

Tanita Corporation

14-2, 1-Chome, Maeno-Cho
Itabashi-Ku, Tokyo, Japan 174-8630
Tel: +81(0)-3-3968-2123
Fax: +81(0)-3-3967-3766
www.tanita.co.jp
ISO 9001 Certified

Tanita Europe GmbH

Dresdener Strasse 25
D-71065 Sindelfingen
Germany
Tel: +49 (0)-7031-6189-6
Fax: +49 (0)-7031-6189-71
www.tanita.de

Tanita Corporation of America Inc.

2625 South Clearbrook Drive
Arlington Heights, IL 60005 USA
Tel: +1-847-640-9241
Fax: +1-847-640-9261
www.tanita.com

Tanita France S.A.S

Villa Labrouste
68 Boulevard Bourdon
92200 Neuilly-Sur-Seine France
Tel: +33 (0)-1-55-24-99-99
Fax: +33(0)-1-55-24-98-68
www.tanita.fr

Tanita Health Equipment H.K. Ltd.

Unit 301-303 3/F Wing On Plaza
62 Mody Road, Tsimshatsui East
Kowloon, Hong Kong
Tel: +852 2838 7111
Fax: +852 2838 8667

©2007 Tanita Corporation. All Rights Reserved.
BF5227811

GB Instruction Manual

Read this Instruction Manual carefully and keep it for future reference.

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GB Table of Contents

INTRODUCTION	1
PRINCIPLES OF ESTIMATING BODY FAT PERCENTAGE	2
FEATURES AND FUNCTIONS	4
PREPARATIONS BEFORE USE	5
INSTALLING THE DISPLAY BOX	5
HANDLING TIPS	6
USEFUL FEATURES	6
AUTOMATIC SHUT-DOWN FUNCTION	6
CHANGING THE WEIGHT MODE	7
GETTING ACCURATE READINGS	7
HOW TO DETERMINE BODY FAT PERCENTAGE	8
TAKING WEIGHT MEASUREMENT ONLY	9
TROUBLESHOOTING	9

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INTRODUCTION

Note: Read this Instruction Manual carefully and keep it handy for future reference.

Thank you for selecting a Tanita Body Fat Monitor/Scale. This model uses the BIA (Bioelectrical Impedance Analysis) technique, a state-of-the-art technology for body fat assessment.

⚠ Safety Precautions

- **Persons with implanted electronic medical equipment, such as a pacemaker, should not use the Body Fat Monitor feature on this Tanita Body Fat Monitor/Scale. This Body Fat Monitor/Scale passes a low-level electrical signal through the body, which may interfere with the operation of a pacemaker.**
- The Tanita Body Fat Monitor/Scales are intended for home use only. This unit is not intended for professional use in hospitals or other medical facilities; it is not equipped with the quality standards required for heavy usage experienced under professional conditions.
- Do not use this unit on slippery surfaces such as wet floors.

Important Notes for Users

This Body fat monitor is intended for adults and children (ages 7-17) with inactive to moderately active lifestyles and adults with athletic body types.

Tanita defines “athlete” as a person involved in intense physical activity of approximately 10 hours per week and who has a resting heart rate of approximately 60 beats per minute or less. Tanita’s athlete definition includes “lifetime of fitness” individuals who have been fit for years but currently exercise less than 10 hours per week.

The body fat monitor function is not intended for pregnant women, professional athletes or bodybuilders.

Recorded data may be lost if the unit is used incorrectly or is exposed to electrical power surges. Tanita takes no responsibility for any kind of loss caused by the loss of recorded data.

Tanita takes no responsibility for any kind of damage or loss caused by these units, or any kind of claim made by a third person.


Note: Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

PRINCIPLES OF ESTIMATING BODY FAT PERCENTAGE

Tanita Body Fat Monitor/Scales allow you to determine your body fat percentage at home as easily as you measure your weight.

What is Body Fat Percentage?

Body fat percentage is the percentage of fat in your body. Too much body fat has been linked to conditions such as high blood pressure, heart disease, diabetes, cancer, and other disabling conditions.

¹ Based on NIH/WHO BMI Guidelines.

² As reported by Gallagher, et al, at NY Obesity Research Center.

To determine the percentage of body fat that is appropriate for your body, consult your physician.

The BIA Method

Tanita Body Fat Monitor/Scales use the BIA (Bioelectrical Impedance Analysis) technique. In this method, a safe, low-level electrical signal is passed through the body. It is difficult for the signal to flow through fat in the human body, but easy to flow through moisture in the muscle and other body tissues. The difficulty with which a signal flows through a substance is called impedance. So the more resistance, or impedance, the signal encounters, the higher the body fat reading.

Tanita’s patented “foot-pad” design sends a safe, low-level electrical signal through the body to determine its composition.

Body Fat Percentage Fluctuations in a Day

Hydration levels in the body may affect body fat readings. Readings are usually highest in the early waking hours, since the body tends to be dehydrated after a long night’s sleep. For the most accurate reading, a person should take a body fat percentage reading at a consistent time of day under consistent conditions.

Besides this basic cycle of fluctuations in the daily body fat readings, variations may be caused by hydration changes in the body due to eating, drinking, menstruation, illness, exercising, and bathing. Daily body fat readings are unique to each person, and depend upon one’s lifestyle, job and activities.

What is Total Body Water Percentage?


Total Body Water Percentage is the total amount of fluid in a person's body expressed as a percentage of their total weight. Approximately 50–65% of the weight of a healthy person is water. Water plays a vital role in many of the body's processes and is found in every cell, tissue and organ. Maintaining a healthy total body water percentage will ensure the body functions efficiently and will reduce the risk of developing associated health problems.

How To Measure Total Body Water?

Deuterium Dilution Method is currently used as a typical method to measure a total body water percentage. This instrument estimates the total body percentage by the analysis (Multiple Regression Analysis) based on the parameter found by the BIA method.

A mutuality over 0.8 is obtained between this instrument and the "deuterium dilution method".

Note: Deuterium Dilution Method is in terms of total body water based on the analysis of the dilution of orally ingested deuterium oxide in urine.

Healthy Total Body Water Percentage (%TBW)

The average total body water percentage ranges for a healthy adult are:

Female: 45 to 60%

Male : 50 to 65%

Note: The total body water percentage of a person with a high percentage of body fat may fall below the average healthy ranges.

Daily Rhythm of Total Body Water

Try to measure total body water percentage at the same time of day and under the same conditions. This increases the accuracy of the readings. The best time is in the early evening, before the evening meal, when hydration levels are more stable. Total body water percentage tends to be underestimated to a degree if measured immediately after getting out of bed, not only because of dehydration but also because the fluid in the body will be stored centrally in the body's trunk. As the day progresses this fluid becomes more evenly distributed due to the effect of gravity.

On the contrary, it may be overestimated after exercise or taking a bath because of the movement of body water or a rise in the body temperature despite actual dehydration.


Important Notes

- Take your measurements at the same time of day and under the same conditions.
- Measured results depend on the movement of water in the body and/or changes in body temperature. Factors affecting total body water levels include sleep, physical activity, eating and drinking, illness and medications, bathing patterns and urination.
- Take fluid immediately if you have any symptoms of dehydration, particularly after physical activity or bathing. Symptoms include a dry mouth, headache, nausea, dizziness, skin dryness, a rise in temperature, loss of concentration, etc. If symptoms persist, consult your doctor immediately.
- Elderly people are at increased risk of dehydration as they have less sensitivity to thirst and other rational symptoms.
- This instrument cannot be used to specifically determine an individual's recommended total body water percentage. Look for the long-term changes in total body water percentage and try to maintain a consistent, healthy total body water percentage.


FEATURES AND FUNCTIONS

Measuring Platform

1. Display Screen
2. Ball-of-Foot Electrodes
3. Heel Electrodes


4. Display Screen
5. Weight/Off
6. Arrow (Select) Button
7. On/Set Button
8. Kg/lb Button

Accessories

1. Wall-Mounting Bracket
2. Screws for Bracket
3. AA-Size Batteries (4)


1.


2.

3.

PREPARATIONS BEFORE USE

▲ Inserting the Batteries

Open the battery cover on the back of the measuring platform. Insert the supplied AA-batteries as indicated.

Note: Be sure that the polarity of the batteries is set properly. If the batteries are incorrectly positioned the fluid may leak and damage floors. If you do not intend to use this unit for a long period of time, it is advisable to remove the batteries before storage.

Please note that the included batteries from factory may have decreased energy levels over time.


▲ Positioning the Monitor

Place the measuring platform on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

Note: To avoid possible injury, do not step on the edge of the platform.


INSTALLING THE DISPLAY BOX


The display box should be set on a table or mounted on a wall. Take care that the display does not fall off the table.

Note: Avoid holding the display unit in your hands while measuring yourself, since its weight (0.6 lb or 0.25 kg) will be added to your weight.

To attach the display box to a wall, follow the steps below:

1. Using a screwdriver, attach the bracket to a wall by inserting the two screws into the holes in the bracket – you may need to pre-drill two holes in the wall.

2. Attach the display box to the bracket by inserting the tabs on the bracket into the holes on the back of the display box.

3. Be sure that the display box is fixed firmly on the bracket. Unless securely attached, the display box may fall off.

HANDLING TIPS

This monitor is a precision instrument utilizing state-of-the-art technology. To keep the unit in the best condition, follow these instructions carefully:

- Do not attempt to disassemble the measuring platform.
- Store the platform horizontally, and place the display unit in its wall-mounted bracket, or in another area where it will not be dropped or stepped on.
Do not wrap cord around the platform or display during storage. Doing so may cause damage to the cord.
- Avoid excessive impact or vibration to the unit.
- Place the unit in an area free from direct sunlight, heating equipment, high humidity, or extreme temperature change.
- Never submerge in water. Use alcohol to clean the electrodes and glass cleaner (applied to a cloth first) to keep them shiny; avoid soaps.
- Do not step on the platform when wet.
- Do not drop any objects onto the platform.

USEFUL FEATURES

- Tanita Body Fat Monitor/Scales allow you to measure both weight and body fat percentage simultaneously and easily, simply by stepping on a scale.
- Your Body Fat Monitor/Scale can be set for Adult or Children (ages 7-17), Male or Female.
- In addition, the BF-522W can be set to Adult Athlete mode.
-  button provides quick weight readings.
- The unit emits a helpful “beep” when activated, and also at various stages in the programming and measuring process. Listen for these sounds which will prompt you to step on the unit, look at your readings, or confirm a setting.

AUTOMATIC SHUT-DOWN FUNCTION


The automatic shut-down function shuts off the power automatically in the following cases:

- If you interrupt the measurement process. The power will shut down automatically within 10 to 20 seconds, depending upon the type of operation.
- If an extraordinary weight is applied to the platform.
- During programming, if you do not touch any of the keys or buttons within 60 seconds.*
- After you have completed the measuring process.

***Note:** If the power shuts off automatically, repeat the steps from “HOW TO DETERMINE BODY FAT PERCENTAGE” (p.7).

CHANGING THE WEIGHT MODE

1. The BF-522W will measure weight in kg (Kilograms), lb (Pounds), or st-lb (British Stone Pounds). The weight mode may be changed anytime when the scale is off. Press the **kg/lb** button as shown in the figure. Press this button until your desired weight mode has been selected.

kg(kilograms)

lb(pounds)

st-lb(stone pounds)


2. After making your selection, the weight mode will be saved and the scale will shut down automatically within 5 seconds.

Note : Pressing the kg/lb button one time, will display the current weight mode setting.

Note : The weight mode may be changed while the weight and body fat % is being displayed.

GETTING ACCURATE READINGS

Heels centered on electrodes

Toes may overhang measuring platform

Note: An accurate reading will not be possible if the soles of your feet are not clean, or if your knees are bent or you are in a sitting position.

To ensure accuracy, readings should be taken without clothing and under consistent conditions of hydration. If you do not undress, always remove your socks or stockings, and be sure the soles of your feet are clean before stepping on the measuring platform. Be sure that your heels are correctly aligned with the electrodes on the measuring platform. Don't worry if your feet appear too large for the unit - accurate readings can still be obtained if your toes overhang the platform. It is best to take readings at the same time of day. Try to wait about three hours after rising, eating, or hard exercise before taking measurements. While readings taken under other conditions may not have the same absolute values, they are accurate for determining the percentage of change as long as the readings are taken in a consistent manner. To monitor progress, compare weight and body fat percentage taken under the same conditions over a period of time.


HOW TO DETERMINE BODY FAT PERCENTAGE

1. Set Age

Press the SET button to turn on the power. The display defaults to Age 30. Use the Up/Down buttons to scroll through numbers. When you reach your age, press the SET button. The unit will beep once to confirm.

Note : Range of user age is between 7-99 years.

Note : If you make a mistake or want to turn the unit off before you have finished programming it, press **Weight Off** button to force quit.

2. Select Female or Male

Use the Up/Down buttons to scroll through Female (♀) or Male (♂), Female/Athlete (♀/🏃) and Male/Athlete (♂/🏃) settings, then press the Set button. The unit will beep once to confirm.

3. Specify the height

The display defaults to 5'7.0"(170cm). Use the Up/Down buttons to specify Height and then press the SET button. The unit will beep once to confirm.

Note : Range of user is from 3' 4.0" – 7' 3.0" or 100cm – 220cm.

4. Get your Readings

The unit will beep twice and "0.0" will appear on the display screen. Make sure you are barefoot before stepping on scale.

Step on to the platform when "0.0" appears. Your weight will be shown first. Continue to stand on the platform. "00000" will appear on the display and disappear one by one from left to right. The readings will be shown alternately 3 times as follows:

Note: Do not step off until Body Fat % is shown

Weight

Body fat percentage

Total body water percentage
(Applicable age 18-99)

Note : If you step onto the platform before "0.0" appears, the display will show "Error" and you will not obtain a reading. Furthermore, if you do not step onto the measuring platform within about 60 seconds after "0.0" appears, the power will shut-off automatically.

TAKING WEIGHT MEASUREMENT ONLY

Press the  button when the unit is off.

When 0.0 appears on the display screen, step on the platform. Your weight will appear on the display and continue to flash as long as you remain on the platform.

TROUBLESHOOTING

If the following problems occur ...then...

- **A wrong weight format appears followed by kg, lb, or st-lb**

Refer to “Switching the Weight Mode” on p.6.

- **”Lo” appears on the display, or all the data appears and immediately disappears.**

Batteries are low. When this message appears, be sure to replace the batteries immediately, since weak batteries will affect the accuracy of your measurements. Change all the batteries at the same time with new AA-size batteries.

Note: Your settings will not be erased from the memory when you remove the batteries.

- **”Error” appears while measuring.**

Please stand on the platform keeping movement to a minimum.

The unit cannot accurately measure your weight if it detects movement.

- **The Body Fat Percentage measurement does not appear or ”-----“ appears after the weight is measured.**

Make sure socks or stockings are removed, and the soles of your feet are clean and properly aligned with the guides on the measuring platform.

- **”OL” appears while measuring.**

Reading cannot be obtained if the weight capacity is exceeded.

- **”Err C” appears after power on.**

The problem might be caused by the connection between the display box and the platform.

If this problem happens, please contact Customer Service.

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