
## *Segmental Body Composition Monitor*

**Model: BC-601FS**

**EN Instruction Manual**

Read this Instruction Manual carefully and keep it for future reference.

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## EN Table of Contents

Introduction .....	1
Features and functions .....	2
Preparations before use .....	3
Getting accurate readings .....	3
Automatic shut-down function .....	3
Setting (resetting) date and time .....	4
Setting and storing personal data in memory .....	5
Taking weight and body fat readings .....	6
How to see the reading of each segment .....	8
Using the visceral fat button .....	9
Visceral fat rating .....	9
Visceral fat trend .....	9
Recall function .....	10
Programming the guest mode .....	11
Taking weight measurement only .....	11
Changing the weight mode .....	11
Troubleshooting .....	12

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# INTRODUCTION

Thank you for selecting a Tanita Body Composition Monitor. This monitor is one in a wide range of home healthcare products produced by Tanita.

This Instruction Manual will guide you through the setup procedures and outline the button features.

Additional information on healthy living can be found on our website [www.tanita.com](http://www.tanita.com).

Tanita products incorporate the latest clinical research and technological innovations.

All data is monitored by our Tanita Medical Advisory Board to ensure accuracy.

The TANITA Body Composition Analyzer is indicated for use in the measurement of weight and impedance, and the estimation of body mass index (BMI), total body fat percent, total body water percent and weight, muscle mass (skeletal and smooth), bone mass, visceral fat rating with healthy range, Daily Calorie Intake (DCI), metabolic age, using B IA (Bioelectrical Impedance Analysis).

The device is indicated for use for healthy children 7-17 years old and healthy adults with active, moderately active, to inactive lifestyles.

**Note:** Read this Instruction Manual carefully and keep it handy for future reference.

## Safety Precautions

Note that people with an electronic medical implant, such as a pacemaker, should not use a Body Composition Monitor as it passes a low-level electrical signal through the body, which may interfere with its operation.

Pregnant women should use the weight function only. All other function are not intended for pregnant women.

Do not place this monitor on slippery surfaces such as a wet floor.

Please do not use scale immediately after bathing or strenuous exercise. Scale platform will be slippery when wet. Please wipe platform with clean, dry cloth to remove any liquids.

## Important Notes for Users

This Body composition monitor is intended for adults and children (ages 7-17) standard and athletic body types.

**Thanks to major investment in the latest BIA Technology and sports science research, TANITA has upgraded the Athlete Mode function, making the measurements more accurate and suitable for a wider audience.**

**Who should use Athlete Mode?**

**Adults from the age of 18 years who either:**

- **Train or exercise for 12 hours or more a week and have been doing so for at least six months.**
- **Are body builders**
- **Are professional athletes who want to monitor their progress at home**
- **People who have a lifetime history of fitness – who used to do more than 12 hours a week – but do less now.**

Recorded data may be lost if the unit is used incorrectly or is exposed to electrical power surges. Tanita takes no responsibility for any kind of loss caused by the loss of recorded data.



Tanita takes no responsibility for any kind of damage or loss caused by these units, or any kind of claim made by a third person.

**Note:** Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

# FEATURES AND FUNCTIONS


## Measuring Platform


AA-size Batteries


Memory Card  
\*SD Logo is a trademark.

## PREPARATIONS BEFORE USE

Open as shown.

Direction of the batteries.

### ⚠ Inserting the Batteries

Turn the scale over and remove the battery cover by releasing the latch, insert 4 x AA batteries (included) with the correct polarity as indicated.

**Note:** Be sure that the polarity of the batteries is set properly. If the batteries are incorrectly positioned the fluid may leak and damage floors. If you do not intend to use this unit for a long period of time, it is advisable to remove the batteries before storage.

Please note that the included batteries from the factory may have decreased energy levels over time.

### ⚠ Positioning the Monitor

Place the measuring platform on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

**Note:** To avoid possible injury, do not step on the edge of the platform.

## Handling Tips

This monitor is a precision instrument utilizing state-of-the-art technology.

To keep the unit in the best condition, please follow these instructions carefully:

- Do not attempt to disassemble the measuring platform.
- Store the unit horizontally, and place it so that the buttons will not be pressed accidentally.
- Avoid excessive impact or vibration to the unit.
- Place the unit in an area free from direct sunlight, heating equipment, high humidity, or extreme temperature change.
- Never submerge in water. Use alcohol to clean the electrodes and glass cleaner (applied to a cloth first) to keep them shiny; avoid using soaps.
- Do not step on the platform with wet feet.
- Do not drop any objects onto the platform.

## GETTING ACCURATE READINGS

To ensure accuracy, readings should be taken without clothing and under consistent conditions of hydration. If you do not undress, always remove your socks or stockings, and be sure the soles of your feet are clean before stepping onto the measuring platform. Be sure your heels are correctly aligned with the electrodes on the measuring platform. Don't worry if your feet appear too large for the unit, accurate readings can still be obtained if your toes overhang the platform. It is best to take readings at the same time of day. Try to wait about three hours after rising, eating, or exercising before taking measurements. While readings taken under other conditions may not have the same absolute values, they are accurate for determining the percentage of change as long as the readings are taken in a consistent manner. To monitor progress, compare weight and body fat percentage taken under the same conditions over a period of time.

**Note:** An accurate reading will not be possible if the soles of your feet are not clean, or if your knees are bent or you are in a sitting position.

## AUTOMATIC SHUT-DOWN FUNCTION

The automatic shut-down function shuts off the power automatically in the following cases:

- If you interrupt the measurement process. The power will shut down automatically within 20 to 70 seconds, depending upon the type of operation.
- If an extreme weight is applied to the platform.
- During programming, if you do not touch any of the buttons within 60 seconds.\*
- After you have completed the measuring process.



\***Note:** If the power shuts off automatically, repeat the steps from "Setting and Storing Personal Data in Memory" (p.4).

# SETTING (RESETTING) DATE AND TIME



## Initial setting after batteries are installed

Press  to turn on the unit and set the date and time. Use the  /  buttons to scroll through numbers. Press SET to confirm.



### 1. Set Year (Applicable year: 2009 – 2050)



### 2. Set Month



### 3. Set Day

### 4. Set Hour (24 hour format)

### 5. Set Minute

Setting complete.

## Resetting time

Press and hold  (Clock) button for 3 seconds to turn on the unit. Then follow the initial setting procedure step 1 to 5.

### Note:

- If you change the date to a date prior to the current date, all tracking data for all users will be lost.
- Replacing the batteries will cause the loss of date and time settings.
- 60 seconds or more of inactivity will cause automatic shut off during this process. No data will be saved as a result.
- If you make a mistake or wish to turn the unit off before you complete the programming, press the OFF button to turn the unit off. No data will be saved as a result.
- If you press (not press and hold) the  (Clock) button, you can see the current date and time on the display.

# SETTING AND STORING PERSONAL DATA IN MEMORY

You will receive readings only if data has been programmed into one of the personal data memories or the Guest mode.


The monitor can store personal details for up to 4 people.

Press SET to turn on the unit. Use the  /  buttons to scroll through numbers. Press SET to confirm.



## Note:

- If you don't operate the unit for 60 seconds after the unit has been turned on, the unit will turn off automatically.
- **If you make a mistake or want to turn the unit off before you have finished programming it, press the  OFF to force quit.**

### 1. Select a personal key number (Up to 4 users)



### 2. Set your Birthday (Year/Month/Day).



### 3. Select Female or Male



Female () , Male () , Female/Athlete () , and Male/Athlete ()

### 4. Specify the Height

### 5. Specify the Activity Level

Setting complete.

To use the DCI (Daily Calorie Intake) feature properly, you need to set the Activity Level.

Use the  button to specify Activity Level. Press the SET to confirm.

Level 1 = Inactive lifestyles (little or no exercise)

Level 2 = Moderately active lifestyles (occasional, low intensity exercise)

Level 3 = Adult involved in intense physical activity (see p.1. for Tanita's definition of Athlete)

**Note:** If you select the Athlete Mode, it will automatically set the Activity Level to Level 3.

The unit will beep twice and your programmed data will be displayed to confirm.

The power will then shut down automatically.

## Note:

- When you change the birthday and/or gender, the unit will show "CLr". If you select "YES" using  /  buttons , the existing personal data and tracking data will be deleted. If you select "no" using  /  , no new data will be saved (existing data will be unchanged) then the unit will be shut off.
- In order to delete personal data, press SET and choose a number with  /  , then press and hold SET. The display will show CLr. Choose YES or no with  /  and press SET to confirm.

# TAKING WEIGHT AND BODY FAT READINGS

## 1. Select Personal Data Number key and Step On

Using your toe, tap your pre-programmed Personal Key.  
The unit will beep and the display will show the programmed data.

**The unit will beep again and the display will show “0.0”.**


(You can change your height setting by pressing SET while “0.0” is shown. Specify your height with  / .)

Press SET to save the changed height, and then the unit will show “0.0” for ready to measure again.)

**Now pull out the hand-electrodes, hold with arms straight down and step onto the platform.**

**Note:** Remember no shoes or socks.


The unit will beep twice to complete the measurement.

### Note:

- **Please grip the handle-electrode as shown, using a light-firm pressure.**
- Do not move until measurement is completed.
- If you step onto the platform before “0.0” appears, the display will show “Error” and you will not obtain a reading. Furthermore if you do not step onto the platform within about 30 seconds after “0.0” appears, the power is shut off automatically.
- Keep your arms in straight up and down while making sure that you do not touch your legs or any part of your body with your arms or hands.


Make sure all of your fingers contact electrodes.

Ensure that your arms are full extended and your elbows do not touch your body.

Make sure your feet contact all of the electrodes.

Do not bend your knees.


## 2. Obtain Your Readings

After the measurement, the unit will display all the readings. Step off the scale and press **RESULT** to see a desired reading.

### Note:

- The unit will automatically shut off after 40 seconds of displaying result with inactivity.
- If your age reaches 100 years old or more, only the weight result will be displayed.
- For children (age 7-17), the unit will only display the Weight, BMI and the Body fat %.

1) Weight (Applicable age: 7 and over)

2) Body fat % (Applicable age: 7-99)

9) Visceral fat level (Applicable age: 18-99)

3) Muscle mass (Applicable age: 18-99)

8) Total body water % (Applicable age: 18-99)

4) Bone mass (Applicable age: 18-99)

7) Metabolic age (Applicable age: 18-99)

5) BMI (Applicable age: 7-99)

6) DCI (Applicable age: 18-99)

## HOW TO SEE THE READING OF EACH SEGMENT (Applicable age : 18-99)

Press **SEGMENT** to show the segmental body fat %, and press **SEGMENT** again to show each segmental reading.

Press **RESULT** to show the muscle mass reading while the segmental body fat % is displayed, and press **SEGMENT** again to show each segmental reading.

Press **RESULT** again while the segmental reading is displayed, except total body fat reading and total muscle mass reading, to choose the segmental body fat % or muscle mass to be displayed.

Press **RESULT** to show the other readings, except segmental readings, while the total body fat % reading or total muscle mass reading is displayed.

### Note:

- The unit will automatically shut off after 40 seconds of inactivity.

Total Body Fat %

Body Fat % - Right Arm

Body Fat % - Left Arm

Body Fat % - Right Leg

Body Fat % - Left Leg

Body Fat % - Trunk Part

Total Muscle Mass

Muscle Mass - Right Arm

Muscle Mass - Left Arm

Muscle Mass - Right Leg



Muscle Mass - Left Leg

Muscle Mass - Trunk Part

### Level Indicator

## USING THE VISCERAL FAT BUTTON

Using  enables you to see your visceral fat rating whenever your readings are displayed.

After the measurement, step off the platform and press  while your readings are displayed.

## VISCERAL FAT RATING

**This feature indicates the rating of visceral fat in your body.**

Visceral fat is the fat that is in the internal abdominal cavity, surrounding the vital organs in the trunk (abdominal) area. Research shows that even if your weight and body fat remains constant, as you get older the distribution of fat changes and is more likely to shift to the trunk area especially post menopause. Ensuring you, have healthy levels of visceral fat may reduce the risk of certain diseases such as heart disease, high blood pressure, and the onset of type 2 diabetes.

The Tanita Body Composition Monitor will provide you with a visceral fat rating from 1 – 59.

### Rating from 1 – 12

Indicates you have a healthy level of visceral fat. Continue monitoring your rating to ensure that it stays within the healthy range.

### Rating from 13 – 59

Indicates you have an excess level of visceral fat. Consider making changes in your lifestyle possibly through diet changes or increasing exercise.

**Source:** Data from Columbia University (New York) & Tanita Institute (Tokyo)


### Note:

- Even if you have a low body fat rate, you may have a high visceral fat level.
- For medical diagnosis, consult a physician.

## VISCERAL FAT TREND

The arrow (→) indicates if your visceral fat level is on the increase or decrease.

The trend indication is displayed only when you take measurements at least 5 times in 2 weeks. Please be sure to take measurement everyday under same conditions.

## RECALL FUNCTION

Recall function can be used to obtain the following desired results in sequence by pressing  /  ;

- Last reading
- The average of the results from 1 week ago
- The average of the results from 2 weeks ago
- The average of the results from 1 month ago
- The average of the results from 2 months ago
- The average of the results from 3 months ago
- The average of the results from 4 months ago
- The average of the results from 5 months ago
- The average of the results from 6 months ago

### While the current reading is displayed after taking measurement;

- 1) Press  to show the previous readings.
- 2) Press  to obtain the other previous readings (weight, total body fat%, total muscle mass, bone mass, BMI, DCI, metabolic age, total body water and Visceral fat).
- 3) Press  to show the previous segmental body fat %, and press  again to show each segmental reading.
- 4) Press  to show the previous muscle mass reading while the previous segmental body fat % is displayed, and press  again to show each segmental reading.
- 5) Press  again while the segmental reading is displayed, except total body fat reading and total muscle mass reading, to choose the previous segmental body fat % or muscle mass to be displayed.
- 6) Press  to show the other readings, except segmental readings, while the total body fat % reading or total muscle mass reading is displayed.
- 7) Press  to show the previous visceral fat reading during any previous reading displayed.
- 8) Press SET again to go back to the current result display.

### While the scale is off;

- 1) Press a personal key to turn on the scale, and press  to show the previous readings.
- 2) To see other previous readings, please refer to the above steps 2) – 5).

### Note:

- The unit will automatically shut off after 40 seconds of inactivity.
- **For children (age 7-17), the unit will only display the Weight, BMI and the Body fat %.**

## PROGRAMMING THE GUEST MODE


Guest mode allows you to program the unit for a one-time use without resetting a Personal key.

Press "GUEST/ ", the unit will beep and the display will show "0.0" indicating it is ready to measure.

**Now pull out the hand-electrodes, hold with arms straight down and step onto the platform.**

After the measurement, the unit will be ready to input your personal information (Age, Gender, Height and Activity Level).

Your readings will be displayed once you complete entering your personal information.

### Note:


- The unit will automatically shut off after 40 seconds of displaying result with inactivity.
- **For children (age 7-17), the unit will only display the Weight, BMI and the Body fat %.**

## TAKING WEIGHT MEASUREMENT ONLY

### Select Weight-Only Key

Using your toe, tap the Weight-only key. After 2 or 3 seconds, "0.0" will appear on the display and the unit will beep. When "0.0" appears on the display, step onto the platform.

The body weight value will appear. If you remain on the platform, the weight result will be displayed for about 20 seconds, and then the power will automatically turn off. If you step off the platform, the weight will continue to show for about 20 seconds then the unit will shut off automatically.


## CHANGING THE WEIGHT MODE

**Note: For use only with models containing the weight mode switch.**

You can change the weight unit indicator by using the switch on the back of the hand-electrodes (kg/lb/st lb).

# TROUBLESHOOTING

## If the following problems occur ... then...

-  appears on the display, or all the data appears and immediately disappears.**  
 Batteries are low. When this message appears, be sure to replace the batteries immediately, since weak batteries will affect the accuracy of your measurements. Change all the batteries at the same time with new batteries.  
 Note: Your settings and past memory will not be erased from the memory when you remove the batteries.
-  appears while measuring.**  
 Please stand on the platform keeping movement to a minimum.  
 The unit cannot accurately measure your weight if it detects movement.  
 The Body Fat Percentage is more than 75%, readings cannot be obtained from the unit.
- The Body Fat Percentage reading does not appear.**  
 Make sure socks or stockings are removed, and the soles of your feet are clean and properly aligned with the guides on the measuring platform.
-  appears while measuring.**  
 Reading cannot be obtained if the weight capacity is exceeded.

### Specifications

#### Weight Capacity

150kg (330lb) (23st 8lb)

#### Weight Increments

0.1kg (0.2lb) (0.2lb)

#### Body Fat Increments

0.1%

#### Personal Data

4 memories

#### Power Supply

4 AA batteries (included)

#### Power Consumption

60 mA maximum

#### Measuring Current

50 kHz, 100  $\mu$ A

### Federal Communications Commission and Canadian ICES Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and Canadian ICES-003. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio or television technician for help.

### Modifications

The FCC requires the user to be notified that any changes or modifications made to this device that are not expressly approved by Tanita Corporation may void the user's authority to operate the equipment.

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**TANITA Corporation of America Inc.**
**TANITA (Shanghai) Trading Co. Ltd.**
**Representative****TANITA Europe B.V.**
**TANITA India Private Limited**