TANITA®
IRONMAN®

innerScan **DUAL**

Radio Wireless Segmental Body Composition Platform

Model: BC-1500 plus

(EN)

Instruction Manual

Read this Instruction Manual carefully and keep it for future reference.

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IRONMAN® and the "M-DOT" logo are registered trademarks of World Triathlon Corporation. Official Product of the IRONMAN TRIATHLON. Used here by permission.

EN Table of Contents

Introduction	1
Features and Functions	1
For Your Safety	2
Preparations before use	4
Getting accurate readings	4
Pairing the BC-1500 plus (“Platform”) with software and measurement	5
Indicator light information	5
BC-1500 plus readings in Healthy Edge [®] software	6
Body Composition Guide	7
What Is Body Fat Percentage?	8
What Is Total Body Water Percentage?	9
What Is Visceral Fat Rating?	9
What Is Basal Metabolic Rate (BMR)?	10
WHAT IS DAILY CALORIE INTAKE (DCI)?	10
WHAT IS METABOLIC AGE?	10
WHAT IS MUSCLE MASS?	11
What Is Physique Rating?	11
WHAT IS BONE MASS?	12
What Is Body Mass Index (BMI)?	12
Troubleshooting	13

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Introduction

Thank you for purchasing this Tanita Body Composition Monitor. This monitor is one of a wide range of home health care products produced by Tanita. This Instruction Manual will guide you through the setup procedures and outline the key features of this monitor. Additional information on healthy living can be found on our website www.tanita.com Tanita products incorporate the latest clinical research and technological innovations. All data is monitored by the Tanita Medical Advisory Board to ensure accuracy. The TANITA body Composition Monitor is intended for use in the measurement of weight and impedance, (and) total and segmental body fat percentage, body water, total and segmental muscle mass, physique rating, bone mass, visceral fat rating with healthy range, basal metabolic rate(BMR), daily calorie intake(DCI), metabolic age and heart rate.

Note: Read this Instruction Manual carefully and keep it handy for future reference.

⚠ Safety Precautions

Note that people with an electronic medical implant, such as a pacemaker, should not use a Body Composition Monitor as it passes a low-level electrical signal through the body, which may interfere with its operation. Pregnant women should use the weight function only. All other function are not intended for pregnant women. Do not place this monitor on slippery surfaces such as a wet floor. Recorded data may be lost if the monitor is used incorrectly or is exposed to electrical power surges. Tanita takes no responsibility for any kind of loss caused by the loss of recorded data. Tanita takes no responsibility for any kind of damage or loss caused by the monitor, or any kind of claim made by a third person.

Important Notes for Users

This Body composition monitor is intended for adults and children (ages 7-17) with inactive to moderately active lifestyles and adults with athletic body types.

Tanita defines “athlete” as a person involved in intense physical activity of approximately 10 hours per week and who has a resting heart rate of approximately 60 beats per minute or less. Tanita’s athlete definition includes “lifetime of fitness” individuals who have been fit for years but currently exercise less than 10 hours per week. The body composition monitor function is not intended for pregnant women, professional athletes or bodybuilders.

Recorded data may be lost if the unit is used incorrectly or is exposed to electrical power surges. Tanita takes no responsibility for any kind of loss caused by the loss of recorded data.

Tanita takes no responsibility for any kind of damage or loss caused by these units, or any kind of claim made by a third person.

Note: Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

Features and Functions

Front side	Back side
<p>CD-ROM</p>	<p>Warranty card</p>
	<p>Quick Start Guide</p>
	<p>AA Battery-LR6 x 4</p>

For Your Safety

This section explains precautionary measures to be taken to avoid injury to the users of this device and others, and to prevent damage to property. Please familiarise yourself with this information to ensure safe operation of this equipment.

	Warning	Failure to follow instructions highlighted with this mark could result in death or severe injury.
	Caution	Failure to follow instructions highlighted with this mark could result in injury or damage to property.
		This mark indicates actions that are prohibited.
		This mark indicates instructions that must always be followed.

Warning

 Prohibited	This monitor must not be used by people with pacemakers or other medical implants. This monitor passes a weak electrical current through the body which could interfere with and cause the malfunction of electrical medical implants, resulting in serious harm to the user.
	Interpretation of measurements (e.g. evaluation of measurements and formulation of exercise programmes based on the measurements) must be performed by a professional. Implementing weight loss measures and exercise programmes based on self-analysis could be harmful to your health. Always follow the advice of a qualified professional.
	Do not step on the edge of the platform. Do not jump onto, or on the monitor. Doing so may cause you to fall or slip, resulting in serious injury.
	Do not place this monitor on a slippery surface such as a wet floor. Doing so may cause you to fall or slip, resulting in serious injury. It may also cause water to get inside the monitor resulting in damage to internal components.
	Do not drop any objects onto the monitor. Avoid excessive impact or vibration to the monitor. Doing so may damage the glass, resulting in injury.
	Do not use the monitor if you have any kind of wound or inflammation on any part of your body that comes into direct contact with the monitor.

Caution

 Prohibited	Do not stand on the display screen or control panel.
	Do not insert your fingers into gaps or holes. Doing so may result in injury.
	Do not use rechargeable batteries. Do not use old batteries together with new batteries, or a mix of different types of batteries at the same time. Doing so may cause battery fluid to leak or the batteries to become excessively hot and rupture, resulting in damage to the monitor or injury.
	Never allow an infant or small child to carry the monitor. The child may drop the monitor, resulting in injury.

Prohibited

This monitor is intended for home use only.

This monitor is not intended for professional use including hospitals, medical or fitness facilities. It is not equipped with the components required for heavy usage. Using the monitor in a professional location will void the warranty.

Do not use the monitor if you have a metal allergy.

The electrodes on this monitor are made from stainless steel which may cause an allergic reaction.

Always...

Place the monitor on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

Avoid excessive impact or vibration to the monitor.

Doing so may damage the monitor or cause it to malfunction, or may cause the setting to be erased.

Do not use or store this monitor in the following locations.

- Locations that are exposed to direct sunlight, near to heaters, or exposed directly to air from air conditioning units, etc.
- Locations where there are severe temperature changes or a large amount of vibration.
- Locations with high levels of humidity or moisture.

Doing so may damage the monitor or cause it to malfunction.

Do not step on the monitor when your body is wet (immediately after a shower or bath, etc.)

Be sure to thoroughly dry your body and feet before using the monitor. Failure to do so may cause you to slip or fall, resulting in injury. Also, the monitor cannot take accurate measurements if your body or feet are wet.

Keep away from water.

Do not use boiling water, benzene, thinner, etc. to clean the monitor.

If the monitor becomes dirty, soak a soft cloth in water or neutral household detergent, wring it out thoroughly and wipe the monitor clean, then use a dry cloth to wipe the monitor dry.

Do not attempt to disassemble the monitor. Do not modify this equipment in any way.

Do not wash the monitor in water.

Do not submerge this monitor in water.

Doing so may damage the monitor or cause it to malfunction.

Be sure to insert the batteries with the polarities (+, -) in the correct position.


If the batteries are incorrectly inserted, battery fluid may leak and damage the floor. (If you do not intend to use this monitor for a long period of time, it is advisable to remove the batteries before storage.)

Do not use the monitor near other products such as microwave ovens that emit electromagnetic waves.

Doing so may cause the monitor to malfunction, preventing accurate measurements.

Another person should assist persons with disabilities who may not be able to take a measurement alone.

Preparations before use

⚠ Inserting the Batteries

Open the battery cover on the back of the measuring platform. Insert the supplied AA-batteries as indicated.

Note: Be sure that the polarity of the batteries is set properly. If the batteries are incorrectly positioned the fluid may leak and damage floors. If you do not intend to use this unit for a long period of time, it is advisable to remove the batteries before storage.

Please note that the included batteries from the factory may have decreased energy levels over time.

⚠ Positioning the Monitor

Place the measuring platform on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

Note: To avoid possible injury, do not step on the edge of the platform.


Handling Tips

This monitor is a precision instrument utilizing state-of-the-art technology.

To keep the unit in the best condition, follow these instructions carefully:

- Do not attempt to disassemble the measuring platform.
- Store the unit horizontally.
- Avoid excessive impact or vibration to the unit.
- Place the unit in an area free from direct sunlight, heating equipment, high humidity, or extreme temperature change.
- Never submerge in water. Use alcohol to clean the electrodes and glass cleaner (applied to a cloth first) to keep them shiny; avoid soaps.
- Do not step on the platform when wet.
- Do not drop any objects onto the platform.
- Do not use this monitor while operating transmitters (such as a cell phone) as this may affect your results.
- If not in use for prolonged period of time it is recommended that you remove the batteries.

Getting accurate readings

Heels centered on electrodes

Toes may overhang measuring platform

To ensure accuracy, readings should be taken without clothing and under consistent conditions of hydration. If you do not undress, always remove your socks or stockings, and be sure the soles of your feet are clean before stepping on the measuring platform. Be sure that your heels are correctly aligned on the measuring platform. Don't worry if your feet appear too large for the unit-accurate readings can still be obtained if your toes overhang the platform.

It is best to take readings at the same time of day. Try to wait about three hours after rising, eating, or hard exercise before taking measurements.

When Heart Rate (HR) measurement is available, the monitor starts to measure HR after Body Composition measurement.

HR value should be higher after exercise or taking a bath.

While readings taken under other conditions may not have the same absolute values, they are accurate for determining the percentage of change as long as the readings are taken in a consistent manner. To monitor progress, compare weight and body fat percentage taken under the same conditions over a period of time.


Note: An accurate reading will not be possible if the soles of your feet are not clean, or if your knees are bent or you are in a sitting position.

Pairing the BC-1500 plus (“Platform”) with software and measurement

The platform needs to pair with the compatible software before measurements are attempted.

Compatible Software: “Healthy Edge”


- 1) Before proceeding,
 - Insert the batteries in the platform.
 - Install the software on your PC (The ANT Driver is installed on your PC automatically).
 - Plug the ANT+ stick into the USB port of your PC.
 - Place the platform within 17feet (5m) of your PC.
- 2) Open the software and select your user profile (which is entered before pairing).
- 3) Click “File” from top left menu, and click “Healthy Edge Options”.
- 4) Select “BC-1500 plus” and click “Connect”. Input the 5 digit serial number found on the bottom of the platform, and then click “OK”. The software memorizes the latest inputted serial number of the platform even if you close the software or PC.
- 5) The platform will beep once and the Green Indicator Light on the platform will blink.
- 6) **Step onto the platform and grip the handle-electrodes while the Green Indicator Light is blinking.** The platform will beep twice and the Green light turns solid green after pairing is complete. Your platform is now paired with your PC.
- 7) For the measurement, click “New Measurement” to use the platform directly from the Software.

Note:

- Please grip the handle-electrode as shown, using a light-firm pressure.
- Place your hands on your side, and do not move until measurement is completed.
- Keep your arms straight up and down while making sure that you do not touch your legs or any part of your body with your arms or hands.
- If you step onto the platform before “Blinking Green Light” appears, the software will show error message and you will not obtain a reading. Furthermore if you do not step onto the platform within about 30 seconds after “Blinking Green Light” appears, the platform will shut off automatically.
- Do not pull the hand-electrodes cable past the red mark.
- Do not carry the unit by holding the hand-electrodes.
- Do not release the hand-electrodes until the cable is completely retracted into the unit.

Indicator light information

BC-1500 plus readings in Healthy Edge[®] software

	Adult (18-99 years old)		Child (7-17 years old)
	Standard	Athlete	
Weight	✓	✓	✓
Total and segmental Body Fat %	✓	✓	✓
Body Water %	✓	✓	n/a
Total and segmental Muscle Mass	✓	✓	n/a
Physique Rating	✓	✓	n/a
BMR	✓	✓	n/a
DCI	✓	✓	n/a
Metabolic Age	✓	✓	n/a
Bone Mass	✓	✓	n/a
Visceral Fat Rating	✓	✓	n/a
Heart Rate	✓	✓	✓ (Reference Data)

Warning

Do not use the body composition reading feature of this product if you have a pacemaker or other electronic implanted medical device. The body composition monitor function is not intended for pregnant women, professional athlete or bodybuilders.

Body Composition Guide

Who Can Use a Body Composition Monitor?

This Body Composition Monitor is intended for adults aged 18-99 years. Children aged 5-17 years can use the monitor for weight, body fat percentage, muscle mass and bone mass: the other features are not applicable to children.

Thanks to major investment in the latest BIA Technology and sports science research, TANITA has upgraded the Athlete Mode function to make measurements more accurate and suitable for a wider range of users.

Who should use Athlete Mode?

Adults aged 18 years and over who either:

- Train or exercise for 12 hours or more a week and have been doing so for at least six months.
- Are body builders.
- Are professional athletes who want to monitor their progress at home.
- Have a lifetime history of fitness and who used to do more than 12 hours a week but do less now.

Pregnant women should only use the weight function.

Other functions are not intended for use when pregnant.

This Body Composition Monitor is intended for home use only. It is not intended for professional use in places such as hospitals or medical or fitness facilities. It is not designed for such heavy usage. Using the monitor in this type of professional environment will invalidate the warranty.

This Body Composition Monitor product provides readings for informational purposes only.


This product is not intended to diagnose or treat any disease or abnormalities. Please consult your physician if you have any questions or concerns related to your health.

How Does a Body Composition Monitor Work?

Tanita Body Composition Monitors calculate your body composition using Dual Frequency Bioelectrical Impedance Analysis (BIA). Safe, low-level electrical signals are passed through the body via the Tanita foot pads on the monitor platform. The signal can flow easily through fluids in muscles and other body tissue but meets resistance as it passes through body fat, because body fat only contains a small amount of fluid. This resistance is called impedance. The impedance readings are then entered into medically researched mathematical formulas to calculate your body composition.

The Tanita RD-953 incorporates medical grade Dual Frequency BIA technology, allowing you to have the highest body composition accuracy in the comfort of your home.

Research has shown that using two different bioelectrical impedance frequencies, provides essential data of a person's intracellular and extracellular status. This advanced technology allows greater accuracy when calculating body composition measurements.

When Is the Best Time To Use My Body Composition Monitor?

Your body water levels naturally fluctuate throughout the day and night. Any significant changes in body water may affect your body composition readings. For example, the body tends to be dehydrated after a long night sleep so if you take a reading first thing in the morning your weight will be lower and your body fat percentage higher. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may also cause variations in your hydration levels.

To get the most reliable reading it is important to use your Body Composition Monitor at a consistent time of day under consistent conditions. A good time to take measurements is before your evening meal.



What Is Body Fat Percentage? (Applicable age 7-99)

Body fat percentage is the amount of body fat as a proportion of your body weight.

Reducing excess levels of body fat has shown to reduce the risk of certain conditions such as high blood pressure, heart disease, diabetes and cancer. The chart below shows the healthy ranges for body fat.

Body Fat Ranges for Standard Children ^{1,2}

Body Fat Ranges for Standard Adults ^{3,4}

Underfat
 Standard (-)
 Standard (+)
 Overfat
 Obese

Your Body Composition Monitor automatically compares your body fat percentage reading to the Healthy Body Fat Range chart.

Underfat: below the healthy body fat range. Increased risk for health problems.

Standard (-)/Standard (+): within the healthy body fat percentage range for your age/gender.

Overfat: above the healthy range. Increased risk for health problems.

Obese: high above the healthy body fat range. Greatly increased risk of obesity-related health problems.

¹ Susan Jebb et al. New Body Fat Reference Curves for Children. *Obesity Research* 2004; 12:A156-157

² HD McCarthy et al. Body Fat Reference Curves for Children. *Int J Obes* 2006; 30: 598-602

³ Gallagher D et al. Healthy percentage body fat ranges: an approach for developing guidelines based on body mass index. *Am J Clin Nutr* 2000; 72: 694-701

⁴ Based on NIH/WHO BMI Guidelines

What Is Total Body Water Percentage? (Applicable age 18-99)

Total Body Water Percentage is the total amount of fluid in a person's body expressed as a percentage of their total weight. Water plays a vital role in many of the body's processes and is found in every cell, tissue and organ. Maintaining a healthy total body water percentage ensures the body functions efficiently and reduces the risk of developing associated health problems. Your body water levels naturally fluctuate throughout the day and night. Your body tends to be dehydrated after a long night's sleep and there are differences in fluid distribution between day and night. Eating large meals, drinking alcohol, menstruation, illness, exercising and bathing may cause variations in your hydration levels.

Your body water percentage reading should be used as a guide and should not be used to specifically determine your recommended total body water percentage. It is important to look for long-term changes in total body water percentage and maintain a consistent, healthy total body water percentage.

Drinking a large quantity of water in one sitting will not instantly change your body water level. In fact, it will increase your body fat reading due to the additional weight gain. Please monitor all readings over time to track relative changes.

Every individual is different, but as a guide the average total body water percentage ranges for a healthy adult are:

Female: 45 to 60%

Male: 50 to 65%

Source: Tanita Research

Note: The total body water percentage will tend to decrease as the percentage of body fat increases. A person with a high percentage of body fat may fall below the average body water percentage. As you lose body fat, the total body water percentage should gradually move towards the typical range given above.

What Is Visceral Fat Rating? (Applicable age 18-99)

This function indicates the amount of visceral fat in your body.

Visceral fat is the fat that is in the internal abdominal cavity, surrounding the vital organs in the abdominal area. Research shows that even if your weight and body fat remains constant, as you get older the distribution of fat changes and is more likely to shift to the abdominal area especially post menopause. Ensuring you have healthy levels of visceral fat may reduce the risk of certain diseases such as heart disease, high blood pressure, and the onset of type 2 diabetes.

The Tanita Body Composition Monitor provides you with a visceral fat rating from 1 – 59.

Rating 1 - 12.5

Indicates you have a healthy amount of visceral fat. Continue monitoring your rating to ensure that it stays within this healthy range.

Rating 13 - 59

Indicates you have an excess amount of visceral fat. Consider making changes in your lifestyle by changing your diet or exercising more.

Source : Data from Columbia University (New York) & Tanita Institute (Tokyo)

Note: You may have a high visceral fat level even if you have a low body fat rate.

- Consult a physician for an accurate medical diagnosis.

What Is Basal Metabolic Rate (BMR)? (Applicable age 18-99)

Your Basal Metabolic Rate (BMR) is the minimum level of energy your body needs when at rest to function effectively. This includes the functioning of your respiratory and circulatory organs, neural system, liver, kidneys, and other organs. You continue to burn calories even when sleeping.

About 70% of the calories your body uses every day are used by your basal metabolism. In addition, energy is used when doing any kind of activity. The more vigorous the activity, the more calories are burned. This is because skeletal muscle (which accounts for approximately 40% of your body weight) acts as your metabolic engine and uses a large amount of energy. Your basal metabolism is greatly affected by the amount of muscle you have, so increasing your muscle mass helps increase your basal metabolism.

By studying healthy individuals, scientists have found that people's metabolic rate change as they age. Basal metabolism rises as a child grows. After peaking at the age of around 16 or 17, it typically starts to gradually decrease.

Having a higher basal metabolism increases the number of calories used and helps decrease the amount of body fat. A low basal metabolic rate makes it harder to lose body fat and overall weight.

How Does a Tanita Body Composition Monitor Calculate BMR?

The basic method of calculating Basal Metabolic Rate (BMR) is a standard equation that uses weight and age.

Tanita has conducted in-depth clinical research into the relationship of BMR and body composition resulting in a much more accurate and personalized reading for the user based on impedance measurements. This method has been medically validated using indirect calorimetry (measuring breath composition).*

*Reliability of the equation for Basal Metabolic Rate: At 2002 Nutrition Week: Scientific and Clinical Forum and Exposition

Title: International Comparison: Resting Energy Expenditure Prediction Models:

The American Journal of Clinical Nutrition.

What Is Daily Calorie Intake (DCI)? (Applicable age 18-99)

"Daily Calorie Intake (DCI)" is the sum of calories for basal metabolism, daily activity metabolism (activities including daily household chores), and diet-induced thermogenesis (energy used in connection with digestion, absorption, metabolism, and other eating activities). It is an estimate of how many calories you can consume within the next 24 hours to maintain your current weight.

How does a Tanita body composition monitor calculate DCI?

The basic way of calculating Basal Metabolic (BMR) Rate is a standard equation using weight and age. Tanita has conducted in-depth research into the relationship of BMR and body composition giving a much more accurate and personalized reading for the user based on the impedance measurement. This method has been medically validated using indirect calorimetry (measuring the breath composition).*

*Reliability on equation for Basal Metabolic Rate: At 2002 Nutrition Week: A

Scientific and Clinical Forum and Exposition Title: International Comparison: Resting Energy Expenditure Prediction Models: The American Journal of Clinical Nutrition.

DCI = BMR × Activity Level
Activity Level

	1	2	3
Female	1.56	1.64	1.82
Male	1.55	1.78	2.10

Source: The World Health Organization (WHO)

What Is Metabolic AGE? (Applicable age 18-99)

This feature calculates your BMR and indicates the average age associated with that type of metabolism.

If your BMR Age is higher than your actual age, it is an indication that you need to improve your metabolic rate. Increased exercise will build healthy muscle tissue, which will improve your metabolic age.

You will obtain a reading between 12 and 90. Under 12 will be displayed as "12" and over 90 displayed as "90".

What Is Muscle Mass? (Applicable age 18-99)

This feature indicates the weight of muscle in your body. The muscle mass displayed includes the skeletal muscles, smooth muscles (such as cardiac and digestive muscles) and the water contained in these muscles.

Muscles play an important role as they act as an engine in consuming energy. As your muscle mass increases, your energy consumption increases helping you reduce excess body fat levels and lose weight in a healthy way.

What Is Physique Rating? (Applicable age 18-99)


This feature assesses your physique according to the ratio of body fat and muscle mass in your body.

As you become more active and reduce the amount of body fat, your physique rating also changes accordingly. Even if your weight does not change, your muscle mass and body fat levels may be changing to make you healthier and reduce your risk of certain diseases.

Each person should set their own goal of which physique they want to achieve, and follow a diet and fitness program to meet that goal.

Result	Physique Rating	Explanation
1	Hidden Obese	Small Frame Obese This person seems to have a healthy body type based on physical appearance; however, they have a high body fat % with low muscle mass level.
2	Obese	Medium Frame Obese This person has a high body fat percentage, with a moderate muscle mass level.
3	Solidly-built	Large Frame Obese This person has both a high body fat % and a high muscle mass.
4	Under exercised	Low Muscle & Average Body Fat % This person has an average body fat % and a less than average muscle mass level.
5	Standard	Ave. Muscle & Ave. Body Fat % This person has average levels of both body fat and muscle mass.
6	Standard Muscular	High Muscle & Ave. Body Fat % (Athlete) This person has an average body fat % and higher than normal muscle mass level.
7	Thin	Low Muscle & Low Fat This person has both a lower than normal body fat % and muscle mass level.
8	Thin and muscular	Thin and muscular (Athlete) This person has lower than normal body fat % while having adequate muscle mass.
9	Very Muscular	Very Muscular (Athlete) This person has lower than normal body fat % while having above average muscle mass.

Source: Data from Columbia University (New York) & Tanita Institute (Tokyo)

What Is Bone Mass? (Applicable age 18-99)

This feature indicates the amount of bone (bone mineral level, calcium or other minerals) in the body.

Research has shown that exercise and the development of muscle tissue are related to stronger, healthier bones. While bone structure is unlikely to make noticeable changes in a short period, it is important that you develop and maintain healthy bones by having a balanced diet and plenty of exercise. People worried about bone disease should consult their physician. People who suffer from osteoporosis or low bone densities due to advanced age, young age, pregnancy, hormonal treatment or other causes, may not get accurate estimations of their bone mass.

Below is the result of estimated bone masses of persons aged 20 to 40, who are said to have the largest amounts of bone masses, by weight. (Source: Tanita Body Weight Science Institute)

Please use the below charts as a guide to compare your bone mass reading.

Women: Average of estimated bone mass

Weight (lb)		
Less than 110 lb	110 lb - 165 lb	165 lb and up
4.3 lb	5.3 lb	6.5 lb

Weight (kg)		
Less than 50 kg	50 kg - 75 kg	75 kg and up
1,95 kg	2,40 kg	2,95 kg

Men: Average of estimated bone mass

Weight (lb)		
Less than 143 lb	143 lb - 209 lb	209 lb and up
5.9 lb	7.3 lb	8.1 lb

Weight (kg)		
Less than 65 kg	65 kg - 95 kg	95 kg and up
2,66 kg	3,29 kg	3,69 kg

Note:

- Persons described below may obtain varying readings and should take the values given for reference purposes only.
 - Elderly persons
 - Women during or after menopause
 - People receiving hormone therapy
- “Estimated bone mass” is a value estimated statistically based on its correlation with the fat-free amount (tissues other than the fat). “Estimated bone mass” does not give a direct judgment on the hardness or strength of the bones or the risks of bone fractures. If you have concerns over your bones, you are recommended to consult a specialist physician.

What Is Body Mass Index (BMI)? (Applicable age 7-99)

This function indicates the relationship between height and weight.

The World Health Organization (WHO) considers an index of 18.5 to 25 as optimal.

Troubleshooting

- Do not step onto the platform before the platform beeps once and the Green Indicator Light blinks.
- Do not step off until hearing a double beep or the Green Light remains solid (stops blinking), the measurement is taken, and the readings are displayed on the applicable device (e.g. Garmin Watch, Tanita Remote Display, and/or PC).
- If you do not step onto the platform within 30 seconds after Green Indicator Light begins blinking, the power will shut-off automatically.
- **Interference may occur if multiple platforms are present. Please use only one platform in any given area.**
- **Constant Red Indicator Light**
Make sure socks or stockings are removed, and the soles of your feet are clean and properly aligned on the measuring platform. Make sure to step onto the platform only **after** it has beeped once or the Green Indicator Light is blinking. The unit cannot accurately measure your weight if it detects movement. Please stand on the platform keeping movement to a minimum.
If Body Fat Percentage is more than 75%, readings cannot be obtained from the unit.
Reading cannot be obtained if the weight capacity is exceeded. Weight capacity is 330lb (150kg).
- **Blinking Red Indicator Light**
Batteries are low. When the Red Indicator Light is blinking, replace the batteries immediately, since weak batteries will affect the accuracy of your measurements. Change all the batteries at the same time with new batteries.

FCC WARNING

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Specifications BC-1500 plus

Weight Capacity

330lb (150kg)

Weight Increments

0.2lb/0.1kg

Power Supply

DC 6V (AA Battery-LR6 × 4 included)

Power Consumption

100mA Maximum

Temperature Range of Usage

41°F–95°F/5°C–35°C

Federal Communications Commission and Canadian ICES Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio or television technician for help.

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EP 0 545 014
U.S. Patent No. 5,415,176
U.S. Patent No. 6,473,641
U.S. Patent No. 6,477,409
U.S. Patent No. 6,480,736
U.S. Patent No. 6,487,445
U.S. Patent No. 6,532,385
U.S. Reissue Patent No. 37,954
Other patents pending

This product is ANT+ certified and is compatible with ANT+ display devices and applications that support the following data:

To learn more visit http://www.thisisant.com/modules/mod_product-directory.php

Representative TANITA Europe B.V.

Hoogoorddreef 56-E
1101 BE Amsterdam
The Netherlands
Tel: +31-20-560-2970
Fax: +31-20-560-2988
www.tanita.eu

TANITA India Private Limited

A-502, Mittal Commercial,
Off. M. V. Road (Andheri Kurla Road),
Marol, Andheri-East,
Mumbai 400 059
Tel: +91-22-3192-6107
Fax: +91-22-2859-9143
www.tanita.in

TANITA Corporation of America, Inc.

2625 South Clearbrook Drive
Arlington Heights, IL 60005 USA.
Tel: +1-847-640-9241
Fax: +1-847-640-9261
www.tanita.com

TANITA Health Equipment H.K. Ltd.

Unit 301-303 3/F Wing On Plaza,
62 Mody Road, Tsimshatsui East,
Kowloon, Hong Kong
Tel: +852-2834-3917
Fax: +852-2838-8667
www.tanita.asia

Manufacturer TANITA Corporation

1-14-2, Maeno-Cho, Itabashi-ku,
Tokyo, 174-8630 JAPAN
Tel: +81-3-3968-7048
Fax: +81-3-3968-2661
www.tanita.co.jp
ISO 9001 Certified