

BC-1000PLUS

RADIO WIRELESS BODY COMPOSITION PLATFORM



The highest standard in radio wireless full-body measurement

Up to 15 full-body measurements when used with Healthy Edge® Lite software (8 user limit)

Thousands of units sold worldwide

Free Healthy Edge® Lite software

Multiple displays: D1000plus / Windows / iPhone compatibility

Available in Black or White



BC-1000plus

RADIO WIRELESS BODY COMPOSITION PLATFORM

Tanita brings the most sophisticated health monitoring technology to health and fitness enthusiasts everywhere! The new **FDA Cleared** BC-1000plus Radio Wireless Body Composition Platform offers consumers the ability to wirelessly link body measurement data in conjunction with a compatible receiving device, providing in-depth information about your body. Body composition monitoring is an essential tool for precisely measuring your body's inner health. Allowing you to accurately monitor the impact of your diet and fitness program.

Tanita and Garmin Ltd., the global leader in satellite navigation, have collaborated to bring consumers the ability to securely and wirelessly transmit data to Garmin Watches. When using this wireless technology, the readings of weight, body fat and body water are wirelessly displayed on the watch, while nine full body readings are stored in the watch. When in range of a consumer's computer, the information is automatically transmitted to the computer using the wireless USB Stick and captured with Tanita's Healthy Edge® Software.

The BC-1000plus will communicate with your computer using the wireless USB Stick and Tanita's Healthy Edge® Software to wirelessly transmit your readings within seconds after stepping on the platform. The wireless interface protocol enables ultra low power applications and offers all the benefits of digital wireless communications. The BC-1000plus scale platform features no visible display, this dynamic feature-the first of its kind-will allow for private readings and is perfect for anyone who wants to keep their health and fitness levels private.

Measurements Include:


- ✓ Weight
- ✓ Metabolic age
- ✓ Daily caloric intake (DCI) or Basal metabolic rate (BMR)
- ✓ Body fat %
- ✓ Bone mass
- ✓ Body water %
- ✓ Visceral fat rating
- ✓ Muscle mass
- ✓ Physique rating

Tanita Connectivity

The BC-1000plus *communicates* using ANT+™ Radio Wireless (with pairing) to transmit securely to multiple compatible display devices, such as:




**iPhone, iPad
& iPod Touch**



**Tanita Healthy
Edge® Mobile
Free App**



Specifications:

- ✓ Capacity: 440lb (200 kg)
- ✓ Weight increments: 0.2lb (0.1kg)
- ✓ Body Fat Increments: 0.1%
- ✓ Transparent Electrodes
- ✓ Tempered Glass Platform
- ✓ Guest Mode & Athlete Mode
- ✓ Power Supply: 4 x AA (included)
- ✓ Dimensions: 13.4" x 14.4" x 1.7"
- ✓ Weight: 6.6 lb
- ✓ USB ANT+ Stick (included) 

Compatible accessories:

- ✓ Healthy Edge® Lite Software (PC) (FREE)
- ✓ Healthy Edge® Plus Software (PC)
- ✓ MY TANITA.COM (FREE)
- ✓ Tanita Wi-Fi Lite App (FREE)
- ✓ D-1000 Plus Remote Display
- ✓ D-1100 Plus Remote Display
- ✓ ANT+ to Wi-Fi Network Adapter
 - ✓ iPhone
 - ✓ iPad
 - ✓ iPod Touch

**Use with ANT+ to Wi-Fi Adapter
for extended range!**



Sense ⇨ Measure ⇨ Communicate